

# PRODUCT CATALOG

**AMSINO**

Advancing Healthcare Worldwide®





**AMSure® Enteral Feeding Systems..... 1**

PUGGLE® Enteral Feeding Pump..... 2  
Sentinel™ Enteral Feeding Pump..... 3  
ENFit® and Non-ENFit Sets..... 3  
ENFit and Non-ENFit Syringes..... 3  
ENFit and Non-ENFit Irrigation Kits and Trays..... 5

**AMSafe® Medication Delivery Systems ..... 6**

Gravity Administration Sets..... 7  
    10 Drops/mL..... 7  
    15 Drops/mL..... 9  
    60 Drops/mL..... 13  
AMSafe-3® Selectable Drip Chamber  
    10, 15 & 60 drops/mL..... 14  
Secondary Piggyback Sets..... 15  
Y-type Blood Sets..... 15  
IV Extension Sets..... 16  
    Standard Extension Sets..... 16  
    Microbore Extension Sets..... 17  
    Smallbore Extension Sets..... 17  
    Flow Controller Extension Sets..... 17  
    Pressure Infusion Extension Sets..... 18  
Needle-Free Connectors..... 18  
Bag Spike..... 18  
Needle-Free Vial Adapter..... 18  
IV Pump Administration Sets..... 19  
    Sure-Lok® Needle-Free IV Pump Administration Sets..... 19  
    Pre-Pierced Y-Injection IV Pump Administration Sets..... 20

**AMSure Urological Care ..... 21**

Foley Catheter Trays & Collection..... 22  
    Pre-Connected Foley Catheter Urine Drainage Bag Trays..... 22  
    Add-a-Foley Catheter Trays..... 23  
    Pre-Connected Foley Catheter Urine Meter Kits..... 24  
    Add-a-Foley Urine Meter Trays..... 25  
    Foley Insertion Trays..... 25  
    Urinary Drainage Bags..... 26  
    Urine Meter..... 27  
    ACCU-MEASURE® Pediatric Urine Meter..... 27  
Foley Catheters..... 28  
    Hydrophilic Latex Foley Catheters..... 28  
    Silicone-Coated Latex Foley Catheters..... 29  
    100% Silicone Foley Catheters..... 30  
Intermittent Urinary Catheters..... 31  
    Urethral Insertion Trays and Kits..... 31  
    Closed System Urethral Drainage Kits..... 32  
    Self-Catheterization Kits with R-Polished™ Eye Catheters..... 33  
    PVC Intermittent Urethral Catheter with R-Polished Eyes..... 34  
    PVC Intermittent Urethral Catheter..... 35  
    Red Rubber Latex Urethral Catheter..... 35  
Irrigation Sets..... 36  
    Cystoscopy & TUR Irrigation Sets..... 36  
Ancillary & Patient Care Products..... 37  
    Urinary Leg Bags..... 37  
    Extension Tube for Leg Bags..... 37

Pediatric Urine Collector..... 37  
Irrigation Trays and Kits..... 38  
Syringes..... 39  
Specimen Containers..... 40

**Sterile Water & Saline Products..... 41**

Sterile Water for Inhalation, USP..... 42  
Sterile Water Prefilled Humidifiers, USP..... 42  
Nebulizer Kits..... 69  
Saline Products..... 43  
    PreFilled Saline Flush Syringes..... 43  
    Sodium Chloride Inhalation Solution, (0.9%) USP..... 44  
Device and Wound Irrigation..... 44

**Suction & Waste Fluid Management Systems..... 45**

EZE-VAC® Suction Canister System..... 46  
EZE-VAC 2.0 Suction Canister System..... 47  
RECEPTAL Liners..... 48  
RECEPTAL Twist Liners..... 49  
RECEPTAL Outer Shells..... 50  
ACCU-MEASURE® Critical Measurement Devices..... 51  
RECEPTAL Solidifiers..... 51  
AMSure® Suction Tubing with Connectors..... 52  
Suction Yankauers and Kits..... 53  
SteriCadi™ Suction Device Holster System..... 54  
Suction Catheters and Trays..... 55  
Suction Catheter Kits..... 56  
Accessories..... 57  
Reusable Hardware..... 60  
Regulator Mounts..... 62  
iRECEPTAL® Surgical Waste Fluid Management System..... 63

**AMSure Respiratory Therapy Systems..... 65**

Respiratory Masks and Nasal Cannulas..... 66  
Nebulizer Kits..... 69  
Accessories..... 70  
Aerosol Drainage Bag..... 70  
Heat and Moisture Exchangers..... 71  
Tracheostomy Clean and Care Trays..... 71  
Sterile Water..... 72  
Suction Catheters, Trays and Kits..... 73  
Suction Tubing with Connectors..... 75  
Suction Yankauers and Kits..... 76

**AMSure Surgical & Nursing Systems..... 77**

Suction Tubing with Connectors..... 78  
Suction Yankauer and Kits..... 79  
Suction Device Holster and Kits..... 79  
Specimen Containers..... 80  
Vaginal Speculums..... 80  
Disposable Umbilical Cord Clamp..... 80  
Ear/Ulcer Syringes..... 81  
Cleansing Enema Sets..... 81  
Connectors..... 82

**Index of REFS ..... 83**

# AMSure®

## ENTERAL FEEDING SYSTEMS



## Optimizing Nutritional Delivery

AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The systems are designed to provide best patient outcomes and includes:

<b>PUGGLE® Enteral Feeding Pump</b> .....	<b>2</b>
<b>Sentinel™ Enteral Feeding Pump</b> .....	<b>3</b>
<b>ENFit &amp; Non-ENFit Enteral Feeding Sets</b> .....	<b>3</b>
<b>ENFit &amp; Non-ENFit Enteral Syringes</b> .....	<b>4</b>
<b>ENFit &amp; Non-ENFit Irrigation Kits &amp; Trays</b> .....	<b>5</b>

# PUGGLE® Enteral Feeding Pump

The versatile pump is designed for pediatric and adult ambulatory applications. Simplified program settings make feeding on-the-go a snap.

**Puggle is only available the US.**

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours



**Dimensions:** 5.4 x 5.0 x 2.5 inches (13.7 x 12.7 x 6.4 cm)  
**Weight:** 1.43 lbs (649 g)  
**Operating Range:** 1 to 400 mL/h in 1 mL increments  
**Delivered Volume:** 1 – 9999 mL in 1 mL increments  
**Accuracy:** +/- 5%



Scan the code below to view the Puggle setup video



REF	DESCRIPTION	QTY/UNIT
PG001	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
PE90ST	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit® and Transition Connector	30/Case
PE90B5T	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
PE90B12T	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case

## Sentinel™ Enteral Feeding Pump



The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each

## ENFit® and Non-ENFit Enteral Feeding Sets



ENF1200K

AMSure Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel™ enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	ENFit CONNECTOR	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
		Spike Only- No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/ EntraFlo™	
E-1004		•				•	•	•	30/Case
E-1200				•		•	•	•	30/Case
E-1203				•	•				30/Case
ENF1004K ©	•	•				•	•	•	30/Case
ENF1200K ©	•			•		•	•	•	30/Case
ENF1203K ©	•			•	•				30/Case
ENF500K ©	•		•			•	•	•	30/Case

EntraFlo™ is a registered trademark of Medline Industries, Inc.  
 Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.  
 Sentinel™ is a registered trademark of Alcor Scientific Inc. ENFit® is a registered trademark of GEDSA.

# ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

## Bulb Syringe

- Flexible soft bulb promotes one-handed administration

## Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage

## Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010 ©			•						•		50/Case
AS011 ©			•				•			•	50/Case
AS115 ©		•			•	•				•	50/Case
AS116		•			•			•			30/Case
ENS015 ©	•			•		•				•	50/Case
ENS015NS	•			•				•			30/Case
ENS115 ©	•				•	•				•	50/Case
ENS115NS	•				•			•			30/Case
ENS115T*	•				•			•			30/Case

\*Includes Transition Adapter

# ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



## Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	30/Case
AS127		•		•		•	•	30/Case
ENK127	•	•		•		•	•	30/Case

## Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	20/Case
ENK136	•	•	•		•	•	•	•	20/Case



# AMSafe®

## MEDICATION DELIVERY SYSTEMS



## Safety First

The AMSafe IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:

Gravity Administration IV Sets.....	7
10 Drops/mL.....	7
15 Drops/mL.....	9
60 Drops/mL.....	13
AMSafe-3® Selectable Drip Chamber, 10, 15 & 60 drops/mL.....	14
Secondary Piggyback Sets .....	15
Y-type Blood Sets .....	15
IV Extension Sets .....	16
Standard Extension Sets.....	16
Microbore Extension Sets.....	17
Smallbore Extension Sets.....	17
Flow Controller Extension Sets.....	17
Pressure Infusion Extension Sets .....	18
Needle-Free Connectors .....	18
Bag Spike.....	18
Needle-Free Vial Adapter .....	18
IV Pump Administration Sets .....	19
Sure-Lok® Needle-Free IV Pump Administration Sets.....	19
Pre-Pierced Y-Injection IV Pump Administration Sets.....	20

# Gravity Administration IV Sets

Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex

## 10 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108303	10	83 in. (211 cm)	Sure-Lok® Needle-free Y-site	16 mL	50/Case
					
108304	10	89 in. (226 cm)	Non-vented, Roller Clamp, (2) Sure-Lok Needle-free Y-sites with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
109602	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
					
108306	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
108310	10	83 in. (211 cm)	Vented/Non-vented Spike, Roller Clamp, (2) Sure-Lok Needle-free Y-sites	16 mL	50/Case
					
MRSA10-84	10	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	10.5 mL	50/Case
					
AC3110	10	105 in. (266 cm)	(2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					
MRA10E-95	10	95 in. (241 cm)	Vented/Non-vented, Pre-pierced Y-site, Slide Clamp, Pre-attached Ext Set, Slide Clamp, MicroClave Clear Needle-free Y-site	18.5 mL	50/Case
					

## 10 Drops/mL

### General Primary Sets (continued)


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>108305</b>	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case




<b>MRA10E-104</b>	10	104 in. (264 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-sites, Pre-pierced Y-site, with Pre-attached Ext Set, Slide Clamp, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13.4 mL	50/Case
-------------------	----	------------------	---	---------	---------



<b>MRSA10E-107</b>	10	107 in. (272 cm)	Vented/Non-vented, (2) MicroClave Clear Needle-free Y-sites, Pre-attached Ext Set	13.4 mL	50/Case
--------------------	----	------------------	---	---------	---------




<b>108308</b>	10	83 in. (211 cm)	(2) Sure-Lok Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
---------------	----	-----------------	--	---------	---------




## 10 Drops/mL

### Flow Controller Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AFS102</b>	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok Needle-free Y-site	18.5 mL	50/Case



<b>AFS208</b> ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) AMSafe® Needle-free Y-sites	19.2 mL	50/Case
-----------------	----	------------------	--	---------	---------



<b>AFS208B</b>	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe Needle-free Y-sites, Slide Clamp	19 mL	50/Case
----------------	----	-----------------	---	-------	---------








# Gravity Administration IV Sets

## 10 Drops/mL



### General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>1010001B</b>	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case
					
<b>108302</b>	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
					
<b>MRA102-84</b>	10	84 in. (213 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
<b>MRSA103-94</b>	10	94 in. (239 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (3) MicroClave Clear Needle-free Y-sites	17.5 mL	50/Case
					
<b>MRSA10E-110</b>	10	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, (2) MicroClave Clear with Pre-attached Ext Set	21.1 mL	50/Case
					

## 15 Drops/mL

### General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AA3101</b>	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
<b>158309</b>	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok® Needle-free Y-site	17 mL	50/Case
					

## 15 Drops/mL General Primary Sets








REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MRSA15-100</b>	15	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site	12.4 mL	50/Case
					
<b>158306</b>	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
<b>MRSA15-91</b>	15	91 in. (231 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	11.4 mL	50/Case
					
<b>159602</b>	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
					
<b>MRA15-100</b>	15	100 in. (254 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13 mL	50/Case
					
<b>1510501</b>	15	105 in. (266 cm)	Vented/Non-vented, (2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					
<b>158305</b>	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
<b>MRSA15E-88</b>	15	88 in. (224 cm)	Vented/Non-vented, Pinch Clamp, Pre-pierced Y-site with Pre-attached Ext Set, Pinch Clamp, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
					
<b>158307</b>	15	89 in. (226 cm)	Sure-Lok Needle-free Y-site, Pre-pierced Y-site with Preattached Ext Set, Sure-Lok Needle-free Y-site, Slide Clamp, Pre-pierced Y-site	18.7 mL	50/Case
					
<b>MRSA15E-120</b>	15	120 in. (305 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site, Rotating Male Luer Lock, Drop-In Ext Set MicroClave Clear Y-site	15.8 mL	50/Case
					

# Gravity Administration IV Sets

## 15 Drops/mL









### Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MRSA15E-96</b>	15	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set, Slide Clamp	12.1 mL	50/Case
					
<b>MRA15B-105</b>	15	105 in. (267 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	13.2 mL	50/Case
					
<b>AS3100</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
<b>AS3101</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
<b>AS3102</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
<b>AS3110A</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					
<b>AS3110B ©</b>	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Pre-pierced Y-sites, Slide Clamp	22 mL	50/Case
					





## 60 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
					
608304	60	89 in. (226 cm)	Non-vented, Roller Clamp, (2) Sure-Lok® Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
608310 ©	60	83 in. (211 cm)	Vented/Non-vented Spike, Roller Clamp, (2) Sure-Lok® Needle-free Y-sites	16 mL	50/Case
					
MRSA60-84	60	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA60E-88	60	88 in. (224 cm)	Vented/Non-vented, Pre-pierced Y-site with Pre-attached Ext Set, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
					
MRSA60-100	60	100 in. (254 cm)	Vented/Non-vented, MicroClave Clear Needle-free Y-site, Pre-pierced Y-site	12.5 mL	50/Case
					

## 60 Drops/mL Piggyback Sets with Back Check Valve


- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRA60B-92	60	92 in. (234 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites, Slide Clamp	11.4 mL	50/Case
					
MRSA602-100	60	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	12.5 mL	50/Case
					




## 60 Drops/mL Burette Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRS60B-106	60	106 in. (269 cm)	150 mL Burette Set (w/MicroClave®, Shut-Off), (2) MicroClave Clear Needle-free Y-sites	12.7 mL	25/Case



MRS60B-117	60	117 in. (297 cm)	24-hour 150 mL Burette Set (w/15 Micron Filter, Shut-off) with (3) MicroClave Clear Needle-free Y-sites	13.5 mL	25/Case
------------	----	------------------	--	---------	---------



## AMSafe-3® Selectable Drip Chamber

### 10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
D38301 ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site	27.7 mL	50/Case



D38302 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.2 mL	50/Case
----------	--------------	-----------------	---	---------	---------




D38305 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Back Check Valve with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.5 mL	50/Case
----------	--------------	-----------------	--	---------	---------




# Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AS3515</b> ©	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
<b>AS3500</b> ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case




# Y-Type Blood Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MRSB-94</b>	10	94 in. (239 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Slide Clamp, Roller Clamp, MicroClave® Y-site	57.2 mL	25/Case



<b>MRSB-102</b>	20	102 in. (259 cm)	200 Micron Filter, Roller Clamps, Hand Pump, MicroClave Y-site	20.2 mL	25/Case
-----------------	----	------------------	--	---------	---------



<b>MRSB-105</b>	10	105 in. (267 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Pre-pierced Y-site, Slide Clamp, Roller Clamp, MicroClave Y-site, Pre-attached Ext Set	59.3 mL	25/Case
-----------------	----	------------------	---	---------	---------



# IV Extension Sets

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.



AE3108

## Standard Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
000606	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok® Needle-free Y-site, Slide Clamp	2.5 mL	50/Case
				
MR4T07	6 in. (15 cm)	MicroClave® Clear Needle-free Connector, T-connector	0.86 mL	50/Case
				
AE3107 ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
				
MR4107	7 in. (18 cm)	Removable MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
MR4207	7 in. (18 cm)	Bonded MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
AE3108 ©	8 in. (20 cm)	Removeable Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
				
AE3108B ©	8 in. (20 cm)	Bonded Sure-Lok Needle-free Connector, Pinch Clamp	1.2 mL	100/Case
				
MR4R10	10 in. (25 cm)	Removable MicroClave Clear Needle-free Connector, Pinch Clamp	0.98 mL	50/Case
				
003201	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
				
003202 ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
				
MR4Y32 ©	32 in. (81 cm)	(2) MicroClave Clear Needle-free Connectors, (2) Slide Clamps	5.02 mL	50/Case
				

## IV Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AE0106</b> ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
<b>AE1106</b> ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
<b>AE0107</b> ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok Needle-free Y-site	0.9 mL	50/Case
<b>AE1112</b>	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
<b>AE1138</b> ©	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case

## Smallbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4307</b>	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.24 mL	50/Case
<b>MR4B07</b>	7 in. (18 cm)	Bifuse Ext Set with (2) MicroClave Clear Needle-free Connector, (2) Slide Clamps	0.35 mL	50/Case
<b>MR4TF04</b>	4 in. (10 cm)	Trifuse Ext Set with (3) MicroClave Clear Needle-free Connectors, (3) Slide Clamps	0.4 mL	50/Case

## Flow Controller Extension Sets

Amsino's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AFS107</b> ©	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
<b>AFS109</b>	19 in. (51 cm)	Flow Controller, NeuFlo® Needle-free Injection Site, Slide Clamp	2.83 mL	50/Case
<b>MRSFC-19</b>	19 in. (48 cm)	Flow Controller, MicroClave Clear Needle-free Connector, Slide Clamp	3.2 mL	50/Case

## Pressure Infusion Extension Sets

Amsino's single-patient-use pressure infusion IV extension set is made to withstand up to 400 psig for rapid administration of medication or fluids in an emergency. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4P07</b>	7 in. (18 cm)	Pressure Infusion (400psig) Ext Set with MicroClave® Clear Needle-free Connector, Purple Slide Clamp	0.43 mL	50/Case




## Needle-Free Connectors

Our AMSafe® needle free valves are manufactured to provide both clinicians and patients save procedural outcomes.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4001</b>	MicroClave Clear Needle-free Pressure Infusion (400psig) Connector	0.04 mL	100/Case
<b>AY0200</b> ©	Sure-Lok® Needle-free Connector	0.33 mL	100/Case
<b>AY0206</b>	NeuFlo® Neutral Pressure Valve	0.08 mL	100/Case



## Bag Spike

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MBSS001</b>	MicroClave Needle-free Connector IV Bag Spike	0.34 mL	50/Case



## Needle-Free Vial Adapter









REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MB4V01</b>	MicroClave Clear Needle-free Connector Single Vial Adapter	0.16 mL	50/Case



# IV Pump Administration Sets

## Sure-Lok® Needle-Free IV Pump Administration Sets


- FDA cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and Pre-Pierced Y-injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP092	10	92 in. (234 cm)	Roller Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP092Y1	10	92 in. (234 cm)	1 Needle-free Valve, Roller Clamp, 1 Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP169BY3	60	109 in. (277 cm)	Needle-free Valve, Roller Clamp, 2 Slide Clamps, Back check Valve, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100FY1	10	100 in. (254 cm)	1 Needle-free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100F12Y1	10	100 in. (254 cm)	Sure-lok Needle-Free Valve, Slide and Roller clamp, 1.2 Micron Filter, Rotating Male Luer Lock with Retractable Collar	22 mL	50/Case
					


## Pre-Pierced Y-Injection IV Pump Administration Sets

- 510(k) cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>FP100BI2</b>	10	100 in. (254 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case



<b>FP109BI3</b>	10	109 in. (277 cm)	Back check Valve, 3 Pre-pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
-----------------	----	------------------	--	-------	---------



<b>FP106FI2</b>	10	106 in. (269 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
-----------------	----	------------------	--	-------	---------



## AMSafe® Prefilled Normal Saline Flush Syringe

- 0.9% Sodium Chloride Injection, USP
- 3 mL, 5 mL, 10 mL and 20 mL
- Sterile Field Ready (SFR) and Sterile Fluid Path packaging configurations



For more details and ordering information, please refer to page 55.



## UROLOGICAL CARE

# Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

- Foley Catheter Trays & Collection ..... 22
  - Pre-Connected Foley Catheter Urine Drainage Bag Trays..... 22
  - Add-a-Foley Catheter Trays ..... 23
  - Pre-Connected Foley Catheter Urine Meter Kits ..... 24
  - Add-a-Foley Urine Meter Trays ..... 25
  - Foley Insertion Trays..... 25
  - Urinary Drainage Bags ..... 26
  - Urine Meter ..... 27
  - ACCU-MEASURE® Pediatric Urine Meter ..... 27
- Foley Catheters ..... 28
  - Hydrophilic Latex Foley Catheters ..... 28
  - Silicone-Coated Latex Foley Catheters ..... 29
  - 100% Silicone Foley Catheters ..... 30
- Intermittent Urinary Catheters ..... 31
  - Urethral Insertion Trays and Kits ..... 31
  - Closed System Urethral Drainage Kits ..... 32
  - Self-Catheterization Kits with R-Polished™ Eye Catheters ..... 33
  - PVC Intermittent Urethral Catheter with R-Polished™ Eyes..... 34
  - PVC Intermittent Urethral Catheter ..... 35
  - Red Rubber Latex Urethral Catheter ..... 35
- Irrigation Sets ..... 36
  - Cystoscopy & TUR Irrigation Sets ..... 36
- Ancillary & Patient Care Products..... 37
  - Urinary Leg Bags ..... 37
  - Extension Tube for Leg Bags ..... 37
  - Pediatric Urine Collector ..... 37
  - Irrigation Trays and Kits ..... 38
  - Syringes ..... 39
  - Specimen Containers ..... 40



# Foley Catheter Trays & Collection



## AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure® Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

**Trays include:** 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S ©	14 Fr.	5-10cc	10/Case
	AS89316S ©	16 Fr.	5-10cc	10/Case
	AS89318S ©	18 Fr.	5-10cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89314 ©	14 Fr.	5-10cc	10/Case
	AS89316 ©	16 Fr.	5-10cc	10/Case
	AS89318 ©	18 Fr.	5-10cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr.	5-10cc	10/Case
	AS89316H	16 Fr.	5-10cc	10/Case
	AS89318H	18 Fr.	5-10cc	10/Case

# Foley Catheter Trays & Collection



AS89300

## AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
<b>AS89300</b>	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10cc	10/Case
<b>AS89330</b> ©	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30cc	10/Case

## AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs). (HCPCS codes: See product chart below)

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag



AS89814S

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	<b>AS89814S</b> ©	14 Fr.	5-10cc	10/Case
	<b>AS89816S</b> ©	16 Fr.	5-10cc	10/Case
	<b>AS89818S</b> ©	18 Fr.	5-10cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	<b>AS89814</b> ©	14 Fr.	5-10cc	10/Case
	<b>AS89816</b> ©	16 Fr.	5-10cc	10/Case
	<b>AS89818</b> ©	18 Fr.	5-10cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	<b>AS89814H</b> ©	14 Fr.	5-10cc	10/Case
	<b>AS89816H</b> ©	16 Fr.	5-10cc	10/Case
	<b>AS89818H</b> ©	18 Fr.	5-10cc	10/Case

# Foley Catheter Trays & Collection

## AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. (HCPCS Code A4354)



AS89800

- **Sterile, not manufactured with natural rubber latex**

Trays include: 400 ml Urine Meter Preconnected with 5cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
<b>AS89800</b>	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10cc	10/Case
<b>AS89830</b>	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30cc	10/Case

## AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)



AS880

- **Deep tray basin designed to prevent contamination**
- **Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection**
- **Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex**

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	<b>AS880</b> ©	10cc	20/Case
	<b>AS890</b> ©	30cc	20/Case
BZK Swab Sticks	<b>AS880K</b> ©	10cc	20/Case
	<b>AS890K</b> ©	30cc	20/Case

## AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube. (HCPCS code A4358)

- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex



AS322 Pear



AS326 Oval



AS32600 Low-Profile

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
<b>AS312</b> ©	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
<b>AS322</b> ©	Pear	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
<b>AS332</b> ©	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
<b>AS302</b> ©	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	20/Case
<b>AS32200</b> ©	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
<b>AS32600*</b> ©	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	20/Case
<b>AS326*</b> ©	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	20/Case

\*100% sterile packed in peel pouch

# Foley Catheter Trays & Collection

## AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. (HCPCS code A4358)



- **Anti-Reflux Chamber helps improve patient safety by reducing urine backflow**
- **Twist valve drains at the bottom of the urine meter for easy drainage**
- **Leak-proof seams ensure fluid containment**
- **Needle-free sampling port reduces the risk of needle stick injuries**
- **Air vents help to maintain consistent urine flow delivering safer patient outcomes**
- **Sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS328</b> ©	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case

## ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)



- **Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag**
- **Needle-free sampling port with luer lock access**
- **Leak-proof seams ensure fluid containment**
- **Meter bottle capacity 150 mL**
- **Adjustable bed sheet clamp for securing drainage bag tubing**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS327100</b>	ACCU-MEASURE® Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case

# Foley Catheters

## AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



### Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
<b>5cc</b> Capacity: 5-10cc	<b>AS41014H</b> ®	14 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41016H</b> ®	16 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41018H</b> ®	18 Fr. × 16 in. (40.6 cm)	10/Box
<b>30cc</b> Capacity: 30-40cc	<b>AS42014H</b>	14 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42016H</b> ®	16 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42018H</b> ®	18 Fr. × 16 in. (40.6 cm)	10/Box

### Coude Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
<b>5cc</b> Capacity: 5-10cc	<b>AS41012HC</b>	12 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41014HC</b>	14 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41016HC</b>	16 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41018HC</b>	18 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41020HC</b>	20 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41022HC</b>	22 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS41024HC</b>	24 Fr. × 16 in. (40.6 cm)	10/Box
<b>30cc</b> Capacity: 30-40cc	<b>AS42016HC</b>	16 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42018HC</b>	18 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42020HC</b>	20 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42022HC</b>	22 Fr. × 16 in. (40.6 cm)	10/Box
	<b>AS42024HC</b>	24 Fr. × 16 in. (40.6 cm)	10/Box

# Foley Catheters

## AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
<b>3cc</b> Capacity: 3-5cc	AS40008 ©	8 Fr. × 10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr. × 10.6 in. (27cm)	10/Box
<b>5cc</b> Capacity: 5-10cc	AS41012 ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014 ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016 ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018 ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020 ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022 ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024 ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026 ©	26 Fr. × 16 in. (40.6cm)	10/Box
	AS41028 ©	28 Fr. × 16 in. (40.6cm)	10/Box
AS41030 ©	30 Fr. × 16 in. (40.6cm)	10/Box	
<b>30cc</b> Capacity: 30-40cc	AS42012 ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014 ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016 ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018 ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020 ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS42022 ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS42024 ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS42026 ©	26 Fr. × 16 in. (40.6cm)	10/Box
	AS42028 ©	28 Fr. × 16 in. (40.6cm)	10/Box
AS42030 ©	30 Fr. × 16 in. (40.6cm)	10/Box	



## AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



### Straight Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
<b>5cc</b> Capacity: 5-10cc	AS41012S ©	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014S ©	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016S ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018S ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020S ©	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022S ©	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024S ©	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026S ©	26 Fr. × 16 in. (40.6cm)	10/Box
<b>30cc</b> Capacity: 30-40cc	AS42016S ©	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018S ©	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020S ©	20 Fr. × 16 in. (40.6cm)	10/Box

### Coude Tips

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
<b>5cc</b> Capacity: 5-10cc	AS41012SC	12 Fr. × 16 in. (40.6 cm)	10/Box
	AS41014SC	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016SC	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018SC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS41020SC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS41022SC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS41024SC	24 Fr. × 16 in. (40.6 cm)	10/Box
<b>30cc</b> Capacity: 30-40cc	AS42018SC	18 Fr. × 16 in. (40.6 cm)	10/Box
	AS42020SC	20 Fr. × 16 in. (40.6 cm)	10/Box
	AS42022SC	22 Fr. × 16 in. (40.6 cm)	10/Box
	AS42024SC	24 Fr. × 16 in. (40.6 cm)	10/Box

# Intermittent Urinary Catheters



AS870

## AMSure® Urethral Insertion Trays and Kits

AMSure® Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection**

### PVC Catheter Tray

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

### Red Rubber Catheter Trays

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
<b>AS870</b> ©	14 Fr PVC*	20/Case
<b>AS87114</b> ©	14 Fr Red Rubber Latex	20/Case
<b>AS87118</b> ©	18 Fr Red Rubber Latex	20/Case

\*Not manufactured with natural rubber latex.

## AMSure® Closed System Urethral Drainage Kits

AMSure® Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection**
- **Sterile**

### Red Rubber Catheter Kits

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

### PVC Catheter Kits

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS89214R

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	<b>AS89214R</b> ©	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	<b>AS89214</b>	14 Fr. × 16 in. (40.6cm)	PVC Latex- Free	20/Case
BZK Swab Sticks	<b>AS89214RK</b> ©	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	<b>AS89216RK</b> ©	16 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case

# Intermittent Urinary Catheters



AS85010

## AMSure® Self-Catheterization Kits with R-POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Introducer tip reduces the potential for urinary tract infection (UTI)**
- **Collection bag with built-in thumb holes allows for ease of use**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

SWAB STICK	REF	CATHETER SIZE	QTY/UNIT
PVP Swab Sticks	<b>AS85008</b> ©	8 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85010</b> ©	10 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85012</b> ©	12 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85014</b> ©	14 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85016</b>	16 Fr. x 16 in. (40.6cm)	50/Case
BZK Swab Sticks	<b>AS85008K</b> ©	8 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85010K</b> ©	10 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85012K</b> ©	12 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85014K</b> ©	14 Fr. x 16 in. (40.6cm)	50/Case
	<b>AS85016K</b> ©	16 Fr. x 16 in. (40.6cm)	50/Case

## AMSure® PVC Intermittent Urethral Catheter with R-POLISHED™ Eyelets

AMSure PVC Intermittent Catheters with R-Polish technology use radio frequency to create ultra-smooth eyelets for optimal patient comfort. They are designed to provide a smooth surface in a catheter that is firm enough for insertion, yet soft enough to be slightly pliable. These urethral catheters are designed for intermittent or short-term use and are available in various French sizes, as well as female, pediatric and male lengths. A Coudé tip is also available to assist in more difficult catheterizations where strictures or other anomalies may be present in the urethra. (HCPCS code A4351)

- **Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal**
- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Designed to minimize trauma and reduce risk of infection**
- **PVC material eliminates latex allergy concerns**
- **Clear shaft provides easy viewing as you measure**
- **Short-term catheter use helps decrease infections**
- **Individually packaged, sterile**
- **BPA-free, not manufactured with natural rubber latex**



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	<b>AS860608</b> ©	8 Fr, 6 in. (15.24 cm)	50/Case
	<b>AS860610</b> ©	10 Fr, 6 in. (15.24 cm)	50/Case
	<b>AS860612</b> ©	12 Fr, 6 in. (15.24 cm)	50/Case
	<b>AS860614</b> ©	14 Fr, 6 in. (15.24 cm)	50/Case
	<b>AS860616</b> ©	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	<b>AS861006</b> ©	6 Fr, 10 in. (25.40 cm)	50/Case
	<b>AS861008</b> ©	8 Fr, 10 in. (25.40 cm)	50/Case
	<b>AS861010</b> ©	10 Fr, 10 in. (25.40 cm)	50/Case
Male	<b>AS861608</b> ©	8 Fr, 16 in. (40.64 cm)	50/Case
	<b>AS861610</b> ©	10 Fr, 16 in. (40.64 cm)	50/Case
	<b>AS861612</b> ©	12 Fr, 16 in. (40.64 cm)	50/Case
	<b>AS861614</b> ©	14 Fr, 16 in. (40.64 cm)	50/Case
	<b>AS861616</b> ©	16 Fr, 16 in. (40.64 cm)	50/Case
	<b>AS861618</b> ©	18 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	<b>AS861610C</b>	10 Fr, 16" (40.65 cm)	30/Case
	<b>AS861612C</b>	12 Fr, 16" (40.65 cm)	30/Case
	<b>AS861614C</b>	14 Fr, 16" (40.65 cm)	30/Case
	<b>AS861616C</b>	16 Fr, 16" (40.65 cm)	30/Case
	<b>AS861618C</b>	18 Fr, 16" (40.65 cm)	30/Case

# Intermittent Urinary Catheters

## AMSure Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free




REF	DESCRIPTION	QTY/UNIT
AS44008 ©	8 Fr, 11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr, 11 in. (27.94 cm)	100/Case
AS44012 ©	12 Fr, 16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr, 16 in. (40.64 cm)	100/Case
AS44016 ©	16 Fr, 16 in. (40.64 cm)	100/Case
AS44018 ©	18 Fr, 16 in. (40.64 cm)	100/Case
AS44020 ©	20 Fr, 16 in. (40.64 cm)	100/Case
AS44022 ©	22 Fr, 16 in. (40.64 cm)	100/Case
AS44024 ©	24 Fr, 16 in. (40.64 cm)	100/Case
AS44028 ©	28 Fr, 16 in. (40.64 cm)	100/Case

# Irrigation Sets

## AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
<b>ABL400</b> ©	Single Lead	20/Case



REF	TYPE	QTY/UNIT
<b>ABL402</b> ©	Two Leads	20/Case



REF	TYPE	QTY/UNIT
<b>ABL404</b> ©	Four Leads	20/Case

## Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex



AS307N

SIZE	REF	DRAIN PORT STYLE	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER & PROTECTIVE CAP	QTY/UNIT
600 mL	AS306N ©	Push-Pull		48/Case
	AS308N ©	Twist-Turn		48/Case
	LBE306N	Push-Pull	•	48/Case
	LBE308N	Twist-Turn	•	48/Case
900 mL	AS307N ©	Push-Pull		48/Case
	AS309N ©	Twist-Turn		48/Case
	LBE307N	Push-Pull	•	48/Case
	LBE309N	Twist-Turn	•	48/Case

## Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS311 ©	18 in. (45.72 cm) Extension Tube	50/Case

## Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS409 ©	200 mL Pediatric Urine Bag	50/Box



# Irrigation Trays and Kits

AMSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need. (HCPCS code A4320)

## Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector



AS121

## Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap



AS126

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector

## Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

## Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		20/Case



AS136

# Ancillary & Patient Care Products

## Syringes

(HCPCS code A4322)

### Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

### Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

### Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



AS010

AS016NS

AS115

TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS011 ©	Poly Pouch	Sterile (EO)	•		50/Case
	AS010 ©	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	AS015 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016NS	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	AS115 ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116	Resealable Pole Bag	Non-Sterile	•	•	30/Case

## AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples. (HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	PACKAGING	STERILITY	LABELS	QTY/UNIT
<b>AS340</b> ©	Poly Pouch	Interior	Patient Information	100/Case
<b>AS341</b> ©	Bulk	Interior	Tamper Evident	100/Case
<b>AS343</b>	Bulk	Non-Sterile	None	500/Case
<b>AS345</b> ©	Form-Fill & Seal	Fully Sterile	None	100/Case



AS340



# AMSure®

## STERILE WATER AND SALINE PRODUCTS



## Inhalation, Irrigation, and Injection Solutions

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

<b>Sterile Water for Inhalation, USP</b> .....	<b>42</b>
<b>Sterile Water Prefilled Humidifiers, USP</b> .....	<b>42</b>
<b>Nebulizer Kits</b> .....	<b>69</b>
<b>Saline Products</b> .....	<b>43</b>
<b>Device and Wound Irrigation</b> .....	<b>44</b>

# AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

## Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	AS1565E ©	500 mL Sterile Water for Inhalation, USP	12/Case
AS1765	AS1765E ©	750 mL Sterile Water for Inhalation, USP	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP	12/Case
AS1566	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1066	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

## Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	1000/Case

## Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

# Saline Products

## AMSafe® Prefilled Saline Flush Syringe

- 3 mL 5 mL, 10 mL and 20 mL syringes, 0.9% Sodium Chloride Injection, USP
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared
- Sterile, non-pyrogenic, DEHP-free, PVC-free, not manufactured with natural rubber latex



## 12 mL and 20 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
<b>IVF1203</b> ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1205</b> ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1210</b> ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF2020</b> ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

### Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
<b>IVF1203SS</b> ©	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1205SS</b> ©	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1210SS</b> ©	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF2020SS</b> ©	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)



## 10 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
<b>IVF1003</b> ©	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1005</b> ©	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
<b>IVF1010</b> ©	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
<b>IVF1010TM</b>	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

## AMSure® Sterile Sodium Chloride Inhalation Solution (vials)

- 3 mL, 5 mL and 15 mL Sterile Sodium Chloride Inhalation Solution (0.9%), USP
- Single-use vials with graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0039	ASR0039E ©	3 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0059	ASR0059E ©	5 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0159	ASR0159E ©	15 mL Sterile Sodium Chloride Inhalation Solution (0.9%) USP	144/Case



ASR0039



ASR0059



ASR0159

## AMSure Device and Wound Irrigation

- 100 mL sterile water or saline bottle with screw cap and inner seal
- Easy open twist screw top container
- Used for device and wound irrigation

REF	DESCRIPTION	QTY/UNIT
AS399 ©	100 mL Sterile Water for Device and Wound Irrigation	48/Case
AS398 ©	100 mL Sterile Saline 0.9% for Device and Wound Irrigation	48/Case



AS399



AS398



## AMSure Sal-Sense™ Single Use Saline Topical Solution 0.9% Sodium Chloride

- Single patient use avoids cross contamination
- Save time from standard practice of drawing solution
- Twist design to open effortlessly
- Easy to irrigate area in any position



REF	DESCRIPTION	QTY/UNIT
AS0159	15 mL Single Use Saline Topical Solution	144/Case–3 Boxes (48 per Box)
AS0359	35 mL Single Use Saline Topical Solution	120/Case–6 Boxes (20 per Box)
AS0559	50 mL Single Use Saline Topical Solution	120/Case–6 Boxes (20 per Box)

# RECEPTAL®

## SUCTION & WASTE FLUID MANAGEMENT SYSTEMS



## Waste Fluid Management Solutions

Amsino is a leading global provider of suction fluid management systems. The breadth of Amsino's RECEPTAL suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:


<b>EZE-VAC® Suction Canister System</b> .....	46
<b>EZE-VAC 2.0 Suction Canister System</b> .....	47
<b>RECEPTAL Liners</b> .....	48
<b>RECEPTAL Twist Liners</b> .....	49
<b>RECEPTAL Outer Shells</b> .....	50
<b>ACCU-MEASURE® Critical Measurement Devices</b> .....	51
<b>RECEPTAL Solidifiers</b> .....	51
<b>AMSure® Suction Tubing with Connectors</b> .....	52
<b>Suction Yankauers and Kits</b> .....	53
<b>SteriCadi™ Suction Device Holster System</b> .....	54
<b>Suction Catheters and Trays</b> .....	55
<b>Suction Catheter Kits</b> .....	56
<b>Accessories</b> .....	57
<b>Reusable Hardware</b> .....	60
<b>Regulator Mounts</b> .....	62
<b>iRECEPTAL® Surgical Waste Fluid Management System</b> .....	63



# EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01 ©	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01 ©	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01 ©	3300 mL Canister, Lid and 2 Elbows	12/Case

# EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.

- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
43210-01	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID × 6 ft. (5 mm × 1.8 m) + 3/16 in ID × 18 in. (5 mm × 46 cm)	30/Case

# RECEPTAL Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
1000	BULK PACKAGED		
	43010-05 ©	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	50/Case
1500	BULK PACKAGED		
	43025-01 ©	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43015S ©	Liner with VAC-GARD and 45g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	50/Case
0L654-01 ©	Liner with VAC-GARD and Pour Spout	50/Case	
2000	BULK PACKAGED		
	43027-01 ©	Liner with Pour Spout,	50/Case
	43041-01 ©	Liner with Ortho Port	50/Case
	43042-01 ©	Liner with VAC-GARD	50/Case
	43044-05 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43020S ©	Liner with VAC-GARD and 60g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	50/Case
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	0G905-01 ©	Liner with Ortho Port	50/Case
0P507-01 ©	Liner with Pour Spout	50/Case	
3000	BULK PACKAGED		
	0L656-01*	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	0L655-01*	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	50/Case
	43030S* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	50/Case

\* Dual port allows for dual suctioning during procedures.  
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID

# RECEPTAL® Twist Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.

- **Premeasured, preloaded and contact-free Solidifier Dispenser Cap**
- **Integrated solidifier eliminates the need for open ports**
- **All-in-one solidifier and liner solution reduces costs associated with extra inventory**
- **100% leak tested seamless liners minimize risk of contamination**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**



43121S



43120S shown above  
in canister 43445-01  
(canister not included)

**Premeasured, preloaded and contact-free solidifier bottle**

**Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure**

**Designed with permanently attached lids and a multi-functional lid tube which helps decrease the risk of fluid exposure**

**VAC-GARD® option, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system**

REF	DESCRIPTION	QTY/UNIT
<b>43110S</b>	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
<b>43115S</b> ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
<b>43116S</b> ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case
<b>43120S</b> ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
<b>43121S</b> ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case

# RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex




Color-coded by size for easy selection.

REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	10/Case
43445-01 ©	2000 mL Outer Shell	10/Case
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

# ACCU-MEASURE® Critical Measurement Devices

RECEPTAL ACCU-MEASURE Device



ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL® ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
43306-01	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case

## RECEPTAL® Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- No shaking required
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4320SP

### RECEPTAL Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
4305SP	500 Solidifier Pouch, solidifies 500 mL	96/Case
4312SP	1200 Solidifier Pouch, solidifies 1200 mL	96/Case
4315SP	1500 Solidifier Pouch, solidifies 1500 mL	96/Case
4320SP	2000 Solidifier Pouch, solidifies 2000 mL	96/Case
4330SP	3000 Solidifier Pouch, solidifies 3000 mL	96/Case



4315SB

4320SB

### RECEPTAL Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
4305SB	500 Solidifier Bottle, solidifies 500 mL	96/Case
4312SB	1200 Solidifier Bottle, solidifies 1200 mL	96/Case
4315SB	1500 Solidifier Bottle, solidifies 1500 mL	96/Case
4320SB	2000 Solidifier Bottle, solidifies 2000 mL	96/Case
4330SB	3000 Solidifier Bottle, solidifies 3000 mL	96/Case
43160SB	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case

# AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex




AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	<b>AS820</b> (50/Case)	<b>AS820B</b> (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> (50/Case)	<b>AS821B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> (50/Case)	<b>AS822B</b> (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> (20/Case)	<b>AS823B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> (20/Case)	<b>AS824B</b> (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> (50/Case)	<b>AS825B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> (50/Case)	<b>AS826B</b> (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> (20/Case)	<b>AS827B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> (20/Case)	<b>AS828B</b> (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> (50/Case)	<b>AS825BM</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> (50/Case)	<b>AS826BM</b> (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> (20/Case)	<b>AS827BM</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> (20/Case)	<b>AS828BM</b> (50/Case)

# Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



## Yankauers

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
	•		•	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
		•		<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
		•	•	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	•			<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
	•		•	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
		•		<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
		•	•	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

## Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

\*Made to order




# SteriCadi™ Suction Device Holster System

SteriCadi is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
SIH0001 ©	SteriCadi Suction Device Holder, Sterile	30/Case



Flex-Field Angle is designed to house a variety of suction devices to help support storage standardization

Suspended tip stays dry

Holder clamp attaches easily to bed rail

Tubing clasp provides a secure attachment

# Suction Catheters and Trays

## Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle tip port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

AS361



REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS361</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS361C</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS362</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS362C</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS363</b> ©	10 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS363C</b> ©	10 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS364</b> ©	12 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS364C</b> ©	12 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS365</b> ©	14 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS365C</b> ©	14 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS366</b> ©	16 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS366C</b> ©	16 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS367</b> ©	18 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS367C</b> ©	18 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case

## Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



SCT10

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
<b>SCT08</b>	Tray, Vinyl Gloves and 17 ½ in. Suction Catheter	8 Fr	DeLee Tip	100/Case
<b>SCT10</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	10 Fr	DeLee Tip	100/Case
<b>SCT12</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	12 Fr	Whistle Tip	100/Case
<b>SCT14</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	14 Fr	Whistle Tip	100/Case
<b>SCT16</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	16 Fr	Whistle Tip	100/Case

# Suction Catheter Kits

AMSure® Suction Catheter Kits provide hassle-free and convenient extraction of mucus or saliva from the upper respiratory airway.

- Uniquely designed two-eyed whistle tip offers safe and effective suctioning
- Sterile, BPA-free, not manufactured with natural rubber latex



AS371

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
<b>AS371</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	1 Vinyl Glove	50/Case
<b>AS381</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	2 Vinyl Gloves	50/Case
<b>AS372</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	1 Vinyl Glove	50/Case
<b>AS382</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	2 Vinyl Gloves	50/Case
<b>AS373</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	1 Vinyl Glove	50/Case
<b>AS383</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	2 Vinyl Gloves	50/Case
<b>AS374</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	1 Vinyl Glove	50/Case
<b>AS384</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	2 Vinyl Gloves	50/Case
<b>AS375</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	1 Vinyl Glove	50/Case
<b>AS385</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	2 Vinyl Gloves	50/Case
<b>AS376</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	1 Vinyl Glove	50/Case
<b>AS386</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	2 Vinyl Gloves	50/Case
<b>AS387</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	18 Fr	2 Vinyl Gloves	50/Case

# Accessories



## RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
43303-01	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
43303-02	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	24/Case
43303-05	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	24/Case



## External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
43309-01	External Shut-Off Valve, 1/4 in. (0.64 cm) Male x Female Elbow	10/Case



## RECEPTAL® Accu-Measure® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
43305-10	Canister Attachment Bracket for 250 mL Canister	10/Case

## Suction Tubing

For permanent Outer Shell-to-vacuum source, or Outer Shell-to-Outer Shell connections.

REF	DESCRIPTION	QTY/UNIT
<b>43455-01</b>	100 ft. (30.48 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	1/Case
<b>43457-25</b>	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case



## EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.

- Non-sterile, BPA-free, not manufactured with natural rubber latex


REF	DESCRIPTION	QTY/UNIT
<b>43317-25</b>	18 in. (0.45 m) Tandem Tubing	25/Case



## RECEPTAL Outer Shell Connection Tubing

RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum "T" of the RECEPTAL canisters between them for easy tandem connection.

REF	DESCRIPTION	QTY/UNIT
<b>CTM50B</b> ©	1/4 in. x 20 in. (6.35 mm x 500 mm), Female Connector	100/Case
<b>CTM20B</b> ©	1/4 in. x 8 in. (6.35 mm x 200 mm), Female Connector	300/Case



## Tubing Manifold

For connection of 2, 3 or 4 canisters.

REF	DESCRIPTION	QTY/UNIT
<b>TMF001</b>	Tubing Manifold	10/Case



# Accessories

## RECEPTAL® Elbow Adapters and Taper

Use with Patient Port.



REF	DESCRIPTION	QTY/UNIT
<b>43318-01</b>	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
<b>43316-01</b>	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
<b>43494-25</b>	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
<b>43473-25</b>	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case

## RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.



REF	DESCRIPTION	QTY/UNIT
<b>44300-05</b>	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	100/Case
<b>43448-25</b>	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case

## 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS807</b> ©	5-in-1 Connector, Sterile	100/Case
<b>AS807B</b>	5-in-1 Connector, Non-sterile	500/Case

## 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- **DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>UYC001S</b>	6-in-1 Universal Y Connector, Sterile	100/Case
<b>UYC001</b>	6-in-1 Universal Y Connector, Non-sterile	500/Case



# Reusable Hardware


## Roll Stands

REF	DESCRIPTION	QTY/UNIT
43474-01	Large head, stainless steel Adjustable Roll Stand, conductive casters, holds up to 12 RECEPTAL® or 8 EZE-VAC® Canisters, use with 1500 mL, 2000 mL and 3000 mL canisters	1/Case
43475-01	Small head, stainless steel Adjustable Roll Stand, conductive casters, holds multiple RECEPTAL and EZE-VAC Canisters	1/Case
43476-P2	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 1500 mL and 2000 mL canisters	1/Case
43476-P3	Built-in Manifold with ON-OFF Valve Adjustable Roll Stand with 4 polycarbonate canister support, use with 3000 mL canisters	1/Case



## Single and In-Line Floor Stands

REF	DESCRIPTION	QTY/UNIT
43431-01	Single Canister Floor Stand with tip-resistant base, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43432-01	Dual Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case
43447-01	In-Line Quad Floor Stand, stainless steel with rollers and conductive casters, use with 1500 mL and 2000 mL RECEPTAL Canisters	1/Case



# Reusable Hardware

## RECEPTAL® Liner System Tandem Canister Supports



REF	DESCRIPTION	QTY/UNIT
43429-01*	For 1500 mL and 2000 mL canisters	1/Case
43439-01	For 3000 mL canisters	1/Case

## RECEPTAL Liner System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43425-05*	For 1000 mL canister, metal	5/Case
43424-10	For 1000 mL canister support, with mounting plate, metal	5/Case
43426-05	For 1500 mL and 2000 mL canisters, metal	5/Case
43438-10	For 3000 mL canisters, with mounting plate for vacuum regulator, metal	10/Case
43428-05*	For 1500 mL and 2000 mL canisters, with mounting plate for vacuum regulator, metal	5/Case
43436-05	For 3000 mL canisters, metal	5/Case
43426P	For 1000 mL, 1500 mL, and 2000 mL canisters, plastic	20/Case
43436P	For 3000 mL canisters, plastic	10/Case
43424-V	For 1000/1500 mL canister support ring with valve	10/Case
43428-V	For 2000 mL canister support ring with valve	10/Case
43438-V	For 3000 mL canister support ring with valve	10/Case

## EZE-VAC® Rigid System Single Canister Supports



REF	DESCRIPTION	QTY/UNIT
43467-01	For EZE-VAC 2.0 800 mL and 1200 mL canisters, metal, small: ID=4.4 in. (11.7 cm), wall to outer edge = 8.5 in. (21.6 cm)	12/Case
43469-12*	For EZE-VAC 1200 mL, 2000 mL and 3300 mL canisters, metal, large: ID=5.6 in. (14.2 cm), wall to outer edge = 9.4 in. (23.9 cm)	12/Case
43469-V	For EZE-VAC 1200 mL, 2000 mL, and 3300 mL with valve	10/Case

\*Shown on Wall Bracket (REF 43444-05)



## Canister Brackets

REF	DESCRIPTION	QTY/UNIT
43444-05	Wall channel bracket for mounting canister, stainless steel	5/Case
43444P	Wall channel bracket for mounting canister, polycarbonate	50/Case
44303-10	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
44304-05	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
RCB001	Metal Rail Bracket with Universal Holder Clamp	30/Case



## Regulator Mounts

REF	DESCRIPTION	QTY/UNIT
43408-05	Use to attach regulator alone or to other supports. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper with 1-1/8 in. (2.86 cm) threaded shank and 2 locknuts	5/Case
43422-05	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper	5/Case
43401-05	Use to attach regulator alone or to other supports. 1/4 in. (0.64 cm) taper with 1-1/4 in. (3.18 cm) shank	5/Case
43400-05	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/4 in. (0.64 cm), male taper	5/Case
OFV001	Metal on-off valve, 1/8 in. NPT	25/Case



Note: Amsino does not sell Regulators

# iRECEPTAL® Surgical Waste Fluid Management System

The iReceptal Digital Surgical Suction System is an all-in-one, closed suction, container, and disposal system. This advanced, easy-to-use surgical suction solution is completely self-contained to provide for independent collection and disposal of waste fluids, along with a docking station for automatic and closed emptying of the rover. iReceptal is currently available in EMEA and China markets. iReceptal is not currently cleared by the US FDA.

- Constantly fully closed system, from collection to disposal
- Staff is never exposed to potentially infectious fluid waste
- Helps optimize surgical turnaround times
- Manifold provides multiple suction lines at the same time with backflow protection
- Separates fluid from tissue and bone fragments
- Waste fluid discharges directly into facility's drain system
- High volume capacity of 20L
- Volume gradation indications
- Easy to maneuver with large wheels
- Lowers disposal costs by eliminating full suction canisters in the waste stream
- Automatically cleans when docked in under 5 minutes

The system is comprised of:

- Mobile Rover
- Docking Station
- Disposable Suction Manifold
- Single-Use Smoke Filter
- Cleaning Detergent




Mobile Rover



Docking Station

REF	DESCRIPTION	QTY/UNIT
<b>iRR100-01</b> ©	iRECEPTAL 1.0 Rover	1/Each
<b>iRR150-01</b> ©	iRECEPTAL 1.5 Rover with 7cm (2.7") LCD Screen	1/Each
<b>iRD102-01</b> ©	iRECEPTAL Docking Station (220VAC)	1/Each
<b>iRD101-01</b>	iRECEPTAL Docking Station (110VAC)	1/Each
<b>iRM100-01</b> ©	4-Port Manifold (disposable)	1/Each
<b>iRC001-01</b>	Detergent for Docking Station	1/Each
<b>iRC002-01</b>	Detergent for Rover	1/Each



### Manifold Storage

- 4 ports for suction tubing with backflow control



### Intelligent Control Panel

- Fluid read-out (Version 1.5 only)

### Completely Closed System

- Avoid direct touch of waste fluid

### Electromagnetic Couplers

- Allow for easy docking




### 20L Volume

- Suitable for large amounts of fluid

### Direct to Drain Systems

- Hands-free disposal and self-cleaning





## Enhancing Oxygen Delivery

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO<sub>2</sub> levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

Respiratory Masks and Nasal Cannulas.....	66
Nebulizer Kits .....	69
Accessories.....	70
Aerosol Drainage Bag .....	70
Heat and Moisture Exchangers .....	71
Tracheostomy Clean and Care Trays .....	71
Sterile Water .....	72
Suction Catheters, Trays and Kits .....	73
AMSure Suction Tubing with connectors .....	75
Suction Yankauers and Kits .....	76

# Respiratory Masks and Nasal Cannulas



## AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

## Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
AS74010	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

# Respiratory Masks and Nasal Cannulas

## High Flow Nasal Oxygen Cannula

Amsino's AMSure® high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most high flow nasal cannula devices
- Single patient use, not manufactured with natural rubber latex



AS75280

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75280</b> ©	High Flow Nasal Oxygen Cannula	Small, Adult	10/Case
<b>AS75285</b>	High Flow Nasal Oxygen Cannula	Medium, Adult	10/Case
<b>AS75282</b> ©	High Flow Nasal Oxygen Cannula	Large, Adult	10/Case

## Nasal Oxygen Cannula

AMSure Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




AS75080

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75080</b> ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
<b>AS75085</b> ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	50/Case

## Cushion Nasal Oxygen Cannula

AMSure Cushion Nasal Oxygen Cannulas with soft foam ear covers provide low-to-medium flow oxygen delivery. Soft, pliable prongs are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex




REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75181</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs with Cushion, Adult	25/Case

## Soft Nasal Oxygen Cannula

AMSure® Soft Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75180

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75180</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Adult	50/Case
<b>AS75190</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Pediatric	50/Case
<b>AS75199</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Infant	50/Case

## Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75010

REF	DESCRIPTION	QTY/UNIT
<b>AS75010</b>	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
<b>AS75020</b>	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case

## Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex



AS75050

REF	DESCRIPTION	QTY/UNIT
<b>AS75040</b>	Adult Aerosol Mask	50/Case
<b>AS75050</b>	Pediatric Aerosol Mask	50/Case

# Respiratory Masks and Nasal Cannulas



## Tracheostomy Mask

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060	Adult Tracheostomy Mask	50/Case



## Nebulizer Kits

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.

- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

### T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

### Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78015 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece	50/Case
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	50/Case
AS78025 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece, Reusable	50/Case



# Accessories

## Respiratory Adapters

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- **Non-sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS730100</b>	Universal Cuff Adapter	22 mm ID Connector	50/Case
<b>AS730200</b>	Multi Adapter	15 mm ID / 22 mm OD Connector	50/Case



## Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- **Pre-connected self-locking plastic chain secures the system to a stable object**
- **Convenient drain port allows for easy one-handed drainage**
- **Non-sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS75035</b>	750 mL Aerosol Drainage Bag	50/Case



# Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
<b>HME001</b> ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
<b>HME002</b> ©	Heat and Moisture Exchangers, Sterile	50/Case

# Tracheostomy Clean and Care Trays



AS861

AMSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS860</b> ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
<b>AS861</b> ©	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	20/Case

# AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

## Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	AS1565E ©	500 mL Sterile Water for Inhalation, USP	12/Case
AS1765	AS1765E ©	750 mL Sterile Water for Inhalation, USP	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP	12/Case
AS1566	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1066	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

## Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	1000/Case

## Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

# Suction Catheters, Trays and Kits

## Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle tip port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS361

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS361</b> ©	6 Fr, 17 ½ in.	Straight Pack	50/Case
<b>AS361C</b> ©	6 Fr, 17 ½ in.	Coiled Pack	50/Case
<b>AS362</b> ©	8 Fr, 17 ½ in.	Straight Pack	50/Case
<b>AS362C</b> ©	8 Fr, 17 ½ in.	Coiled Pack	50/Case
<b>AS363</b> ©	10 Fr, 23 in.	Straight Pack	50/Case
<b>AS363C</b> ©	10 Fr, 23 in.	Coiled Pack	50/Case
<b>AS364</b> ©	12 Fr, 23 in.	Straight Pack	50/Case
<b>AS364C</b> ©	12 Fr, 23 in.	Coiled Pack	50/Case
<b>AS365</b> ©	14 Fr, 23 in.	Straight Pack	50/Case
<b>AS365C</b> ©	14 Fr, 23 in.	Coiled Pack	50/Case
<b>AS366</b> ©	16 Fr, 23 in.	Straight Pack	50/Case
<b>AS366C</b> ©	16 Fr, 23 in.	Coiled Pack	50/Case
<b>AS367</b> ©	18 Fr, 23 in.	Straight Pack	50/Case
<b>AS367C</b> ©	18 Fr, 23 in.	Coiled Pack	50/Case

## Suction Catheter Kits

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Uniquely designed two-eyed whistle tip offers safe and effective suctioning**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS371

REF	DESCRIPTION	FEATURES	GLOVES	QTY/UNIT
<b>AS371</b> ©	6 Fr, 17 ½ in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS381</b> ©	6 Fr, 17 ½ in.	Whistle Tip	2 Vinyl Gloves	50/Case
<b>AS372</b> ©	8 Fr, 17 ½ in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS382</b> ©	8 Fr, 17 ½ in.	Whistle Tips	2 Vinyl Gloves	50/Case
<b>AS373</b> ©	10 Fr, 23 in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS383</b> ©	10 Fr, 23 in.	Whistle Tip	2 Vinyl Gloves	50/Case
<b>AS374</b> ©	12 Fr, 23 in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS384</b> ©	12 Fr, 23 in.	Whistle Tip	2 Vinyl Gloves	50/Case
<b>AS375</b> ©	14 Fr, 23 in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS385</b> ©	14 Fr, 23 in.	Whistle Tip	2 Vinyl Gloves	50/Case
<b>AS376</b> ©	16 Fr, 23 in.	Whistle Tip	1 Vinyl Glove	50/Case
<b>AS386</b> ©	16 Fr, 23 in.	Whistle Tip	2 Vinyl Gloves	50/Case
<b>AS387</b> ©	18 Fr, 23 in.	Whistle Tip	2 Vinyl Gloves	50/Case

## Mini Suction Catheter Trays

The AMSure® Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**




SCT10

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>SCT08</b>	8 Fr, 17 ½ in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
<b>SCT10</b>	10 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
<b>SCT12</b>	12 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
<b>SCT14</b>	14 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
<b>SCT16</b>	16 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case

## Closed Suction Catheter Trays

24-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

- **Closed system to prevent cross infections in facility**
- **Clean easily with needle free port**
- **Avoid interruptions of oxygen supply due to suction**



VS365TP124

VS365TC124

REF	DESCRIPTION	FLEX CONNECTOR	QTY/UNIT
<b>Endotracheal Kits with MDI Port Adapter</b>			
<b>VS363RC124</b>	10 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS364RC124</b>	12 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS365RC124</b>	14 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS363TP124</b>	10 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
<b>VS364TP124</b>	12 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
<b>VS365TP124</b>	14 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
<b>Tracheostomy Kits with MDI Port Adapter</b>			
<b>VS363TC124</b>	10 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS364TC124</b>	12 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS365TC124</b>	14 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
<b>VS363TTP124</b>	10 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
<b>VS364TTP124</b>	12 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
<b>VS365TTP124</b>	14 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case

# Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	<b>AS820</b> © (50/Case)	<b>AS820B</b> © (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> © (50/Case)	<b>AS821B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> © (50/Case)	<b>AS822B</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> © (20/Case)	<b>AS823B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> © (20/Case)	<b>AS824B</b> © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> © (50/Case)	<b>AS825B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> © (50/Case)	<b>AS826B</b> © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> © (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> © (20/Case)	<b>AS827B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> © (20/Case)	<b>AS828B</b> © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> © (1 Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> © (50/Case)	<b>AS825BM</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> © (50/Case)	<b>AS826BM</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> © (20/Case)	<b>AS827BM</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> © (20/Case)	<b>AS828BM</b> © (50/Case)

# Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



AS835

AS832

## Yankauers

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	•			<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
	•		•	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
		•		<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
		•	•	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	•			<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
	•		•	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
		•		<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
		•	•	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

## Pre-Connected Yankauer & Tubing

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

\*Made to order

# AMSure®

## SURGICAL AND NURSING SOLUTIONS



## Everyday Patient Care Solutions

AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

Suction Tubing .....	78
Suction Yankauers and Kits .....	79
Suction Device Holster System .....	79
Specimen Containers.....	80
Vaginal Speculums.....	80
Disposable Umbilical Cord Clamp.....	80
Ear/Ulcer Syringes .....	81
Cleansing Enema Sets.....	81
Connectors.....	82



# Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	AS820 © (50/Case)	AS820B © (100/Case)
	6 ft. (1.8 m)	AS821 © (50/Case)	AS821B © (50/Case)
	10 ft. (3.0 m)	AS822 © (50/Case)	AS822B © (50/Case)
	12 ft. (3.7 m)	AS823 © (20/Case)	AS823B © (50/Case)
	20 ft. (6.1 m)	AS824 © (20/Case)	AS824B © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	AS825 © (50/Case)	AS825B © (50/Case)
	10 ft. (3.0 m)	AS826 © (50/Case)	AS826B © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	AS826C © (50/Case)	—
	12 ft. (3.7 m)	AS827 © (20/Case)	AS827B © (50/Case)
	20 ft. (6.1 m)	AS828 © (20/Case)	AS828B © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	AS829B-01 © (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	AS825M © (50/Case)	AS825BM © (50/Case)
	10 ft. (3.0 m)	AS826M © (50/Case)	AS826BM © (50/Case)
	12 ft. (3.7 m)	AS827M © (20/Case)	AS827BM © (50/Case)
	20 ft. (6.1 m)	AS828M © (20/Case)	AS828BM © (50/Case)

# Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- **Translucent pathway allows for fluid visualization**
- **Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage**
- **Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex**



AS835

AS832

## Yankauers

STYLE	TIP TYPE			STERILE	NON-STERILE
	REGULAR	BULBOUS	VENTED		
Rigid	●			<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
	●		●	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
		●		<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
		●	●	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	●			<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
	●		●	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
		●		<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
		●	●	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

## Yankauer & Tubing Kits

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	20/Case

\*Made to order

# Suction Device Holster System

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- **Suction devices can be isolated to help reduce cross-contamination and improve safety measures**
- **Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>SIH0001</b> ©	SteriCadi Suction Device Holder with Clip	30/Case



# Specimen Containers

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	QTY/UNIT
AS340 ©		•		•		•	100/Case
AS341 ©			•	•	•		100/Case
AS343 ©			•				500/Case
AS345 ©	•			•			100/Case



AS340

# Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- Round-edge design helps to increase patient comfort
- Single-use disposable
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS032S	Small Vaginal Speculum	100/Case
AS032M	Medium Vaginal Speculum	100/Case
AS032L	Large Vaginal Speculum	100/Case



AS032M

# Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- Safety lock is designed to decrease risk of accidental disconnection
- Grooved teeth help provide an effective grip
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UCC100	Umbilical Cord Clamp, sterile	50/Case





## Ear/Ulcer Syringes

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS00501 ©	1 oz. syringe	Non-sterile, Bulk	300/Case
AS00502 ©	2 oz. syringe	Non-sterile, Bulk	200/Case
AS00503 ©	3 oz. syringe	Non-sterile, Bulk	150/Case
AS00501S ©	1 oz. syringe	Sterile, Single Pack	50/Case
AS00502S ©	2 oz. syringe	Sterile, Single Pack	50/Case
AS00503S ©	3 oz. syringe	Sterile, Single Pack	50/Case



AS330

## Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

### Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

### Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS333

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS330 ©	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case

# Connectors

## 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 9/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS807</b> <sup>Ⓒ</sup>	5-in-1 Connector	Sterile	100/Case
<b>AS807B</b>	5-in-1 Connector	Non-sterile, Bulk	500/Case



AS807

## 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>UYC001S</b>	6-in-1 Universal Y Connector, Sterile	100/Case
<b>UYC001</b>	6-in-1 Universal Y Connector, Non-sterile	500/Case



# Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606	16	43318-01	55	AS011	4, 39	AS375	56, 73
003201	16	43400-05	62	AS015	39	AS376	56, 73
003202	16	43401-05	62	AS0159	44	AS381	56, 73
0G901-01	48	43408-05	62	AS016NS	39	AS382	56, 73
0G902-01	48	43422-05	62	AS032L	80	AS383	56, 73
0G903-01	48	43423-01	50	AS032M	80	AS384	56, 73
0G904-01	48	43424-10	61	AS032S	80	AS385	56, 73
0G905-01	48	43424-V	61	AS0352	42, 72	AS386	56, 73
0L654-01	48	43425-05	61	AS0352E	42, 72	AS387	56, 73
0L655-01	48	43426-05	61	AS0359	44	AS398	44
0L656-01	48	43426P	61	AS0552	42, 72	AS399	44
0P507-01	48	43428-05	61	AS0552E	42, 72	AS409	37
0P508-01	48	43428-V	61	AS0559	44	AS40008	29
1010001B	9	43429-01	61	AS1000	42, 72	AS40010	29
107201	7	43431-01	60	AS1065	42, 72	AS41012	29
108301	7	43432-01	60	AS1065E	42, 72	AS41012HC	28
108302	9	43436-05	61	AS1066	42, 72	AS41012S	30
108303	7	43436P	61	AS115	4, 39	AS41012SC	30
108304	7	43438-10	61	AS116	4, 39	AS41014	29
108305	8	43438-V	61	AS121	5, 38	AS41014H	28
108306	7	43439-01	61	AS126	5, 38	AS41014HC	28
108308	8	43444-05	62	AS127	5, 38	AS41014S	30
108310	7	43444P	62	AS130	5, 38	AS41014SC	30
109602	7	43445-01	50	AS136	5, 38	AS41016	29
1510501	10	43447-01	60	AS1565	42, 72	AS41016H	28
158305	10	43448-25	59	AS1565E	42, 72	AS41016HC	28
158306	10	43449-05	50	AS1566	42, 72	AS41016S	30
158307	10	43455-01	58	AS1765	42, 72	AS41016SC	30
158309	9	43457-25	58	AS1765E	42, 72	As41018	29
159602	10	43467-01	61	AS302	26	AS41018H	28
43010-05	48	43469-12	61	AS306N	37	AS41018HC	28
43010S	48	43469-V	61	AS307N	37	AS41018S	30
43015S	48	43473-25	59	AS308N	37	AS41018SC	30
43020S	48	43474-01	60	AS309N	37	AS41020	29
43025-01	48	43475-01	60	AS3100	11	AS41020HC	28
43027-01	48	43476-P2	60	AS3101	11	AS41020S	30
43030S	48	43476-P3	60	AS3102	11	AS41020SC	30
43041-01	48	43494-25	59	AS3110A	11	AS41022	29
43042-01	48	43503-01	50	AS3110B	11	AS41022HC	28
43044-01	48	44300-05	59	AS312	26	AS41022S	30
43044-05	48	44303-10	62	AS322	26	AS41022SC	30
43055B	51	44304-05	62	AS32200	26	AS41024	29
43055P	51	607201	13	AS326	26	AS41024HC	28
43110S	49	608304	13	AS328	27	AS41024S	30
43115S	49	608305	13	AS32600	26	AS41024SC	30
43116S	49	608306	13	AS327100	27	AS41026	29
43120S	49	608310	13	AS330	81	AS41026S	30
43121S	49	AA3101	9	AS333	81	AS41028	29
43125B	51	ABL400	36	AS340	40, 80	AS41030	29
43125P	51	ABL402	36	AS341	40, 80	AS42012	29
43155B	51	ABL404	36	AS343	40, 80	AS42014	29
43155P	51	AC3110	7	AS345	40, 80	AS42014H	28
43160SB	51	AE0106	17	AS3500	15	AS42016	29
43202-01	47	AE0107	17	AS3515	15	AS42016H	28
43202-820821	47	AE1106	17	AS361	55, 73	AS42016HC	28
43203-01	46	AE1112	17	AS361C	55, 73	AS42016S	30
43206-01	46	AE1138	17	AS362	55, 73	AS42018	29
43205B	51	AE3107	16	AS362C	55, 73	AS42018H	28
43205P	51	AE3108	16	AS363	55, 73	AS42018HC	28
43210-01	47	AE3108B	16	AS363C	55, 73	AS42018S	30
43213-01	46	AFS102	8	AS364	55, 73	AS42018SC	30
43303-01	57	AFS107	17	AS364C	55, 73	AS42020	29
43303-02	57	AFS109	17	AS365	55, 73	AS42020HC	28
43303-05	57	AFS208	8	AS365C	55, 73	AS42020S	30
43305-10	51	AFS208B	8	AS366	55, 73	AS42020SC	30
43305-10	57	AS00501	81	AS366C	55, 73	AS42022	29
43306-01	51	AS00501S	81	AS367	55, 73	AS42022HC	28
43309-01	57	AS00502	81	AS367C	55, 73	AS42022SC	30
43305B	51	AS00502S	81	AS371	56, 73	AS42024	29
43305P	51	AS00503	81	AS372	56, 73	AS42024HC	28
43316-01	59	AS00503S	81	AS373	56, 73	AS42024SC	30
43317-25	58	AS010	4, 39	AS374	56, 73	AS42026	29

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
AS42028	29	AS833	53, 76, 79	AS89818S	24	MR4107	16
AS42030	29	AS833B	53, 76, 79	AS89830	25	MR4207	16
AS44008	35	AS834	53, 76, 79	ASR0035	42, 72	MR4307	17
AS44010	35	AS834B	53, 76, 79	ASR0035E	42, 72	MR4B07	17
AS44012	35	AS835	53, 76, 79	ASR0039	44	MR4P07	18
AS44014	35	AS835B	53, 76, 79	ASR0039E	44	MR4R10	16
AS44016	35	AS836	53, 76, 79	ASR0055	42, 72	MR4T07	16
AS44018	35	AS836B	53, 76, 79	ASR0055E	42, 72	MR4TF04	17
AS44020	35	AS837	53, 76, 79	ASR0059	44	MR4Y32	16
AS44022	35	AS837B	53, 76, 79	ASR0059E	44	MRA102-84	9
AS44024	35	AS85008	33	ASR0159	44	MRA10E-104	8
AS44028	35	AS85008K	33	ASR0159E	44	MRA10E-95	7
AS730100	70	AS85010	33	AY0200	18	MRA15-100	10
AS730200	70	AS85010K	33	AY0206	18	MRA15B-105	11
AS74010	66	AS85012	33	CTM20B	58	MRA60B-92	13
AS75010	68	AS85012K	33	CTM50B	58	MRS60B-106	14
AS75020	68	AS85014	33	D38301	14	MRS60B-117	14
AS75035	70	AS85014K	33	D38302	14	MRSA10-84	7
AS75040	68	AS85016	33	D38305	14	MRSA103-94	9
AS75050	68	AS85016K	33	E-0500M	3	MRSA103-94	9
AS75060	69	AS860	71	E-1004	3	MRSA10E-107	8
AS75080	67	AS860608	34	E-1200	3	MRSA15-100	10
AS75085	67	AS860610	34	E-1203	3	MRSA15-91	10
AS75180	68	AS860612	34	ENF1004K	3	MRSA15E-120	10
AS75181	67	AS860614	34	ENF1200K	3	MRSA15E-88	10
AS75190	68	AS860616	34	ENF1203K	3	MRSA15E-96	11
AS75199	68	AS861	71	ENF500K	3	MRSA60-100	13
AS75280	67	AS861006	34	ENK127	5	MRSA60-84	13
AS75282	67	AS861008	34	ENK136	5	MRSA602-100	13
AS75285	67	AS861010	34	ENS015	4	MRSA60E-88	13
AS78010	69	AS861608	34	ENS015NS	4	MRSB-102	15
AS78015	69	AS861610	34	ENS115	4	MRSB-105	15
AS78020	69	AS861610C	34	ENS115NS	4	MRSB-94	15
AS78025	69	AS861612	34	FP092	19	MRSFC-19	17
AS807	59, 82	AS861612C	34	FP092Y1	19	OFV001	62
AS807B	59, 82	AS861614	34	FP100BI2	20	OM001	66
AS820	52, 75, 78	AS861614C	34	FP100BY2	19	PE90B12T	2
AS820B	52, 75, 78	AS861616	34	FP100FY1	19	PE90B5T	2
AS821	52, 75, 78	AS861616C	34	FP100F12Y1	19	PE90ST	2
AS821B	52, 75, 78	AS861618	34	FP106FI2	20	PG001	2
AS822	52, 75, 78	AS861618C	34	FP106FY2	19	RCB001	62
AS822B	52, 75, 78	AS870	31	FP109BI3	20	S-1000-SI	3
AS823	52, 75, 78	AS87114	31	FP109BY3	19	SCT08	55, 74
AS823B	52, 75, 78	AS87118	31	FP169BY3	19	SCT10	55, 74
AS824	52, 75, 78	AS880	25	HME001	71	SCT12	55, 74
AS824B	52, 75, 78	AS880K	25	HME002	71	SCT14	55, 74
AS825	52, 75, 78	AS890	25	iRC001-01	63	SCT16	55, 74
AS825B	52, 75, 78	AS890K	25	iRD101-01	63	SIH0001	54, 79
AS825BM	52, 75, 78	AS89214	32	iRD102-01	63	TMF001	58
AS825M	52, 75, 78	AS89214R	32	iRM100-01	63	UCC100	80
AS825P	53, 76, 79	AS89214RK	32	iRR100-01	63	UYC001	59, 82
AS826	52, 75, 78	AS89216RK	32	iRR150-01	63	UYC001S	59, 82
AS826B	52, 75, 78	AS89300	23	IVF1003	43	VS363RC124	74
AS826BM	52, 75, 78	AS89314	22	IVF1005	43	VS363TC124	74
AS826C	52, 75, 78	AS89314H	22	IVF1010	43	VS363TP124	74
AS826M	52, 75, 78	AS89314S	22	IVF1010TM	43	VS363TTP124	74
AS827	52, 75, 78	AS89316	22	IVF1203	43	VS364RC124	74
AS827B	52, 75, 78	AS89316H	22	IVF1203SS	43	VS364TC124	74
AS827BM	52, 75, 78	AS89316S	22	IVF1205	43	VS364TP124	74
AS827M	52, 75, 78	AS89318	22	IVF1205SS	43	VS364TTP124	74
AS827P	53, 76, 79	AS89318H	22	IVF1205SS	43	VS365RC124	74
AS828	52, 75, 78	AS89318S	22	IVF1210	43	VS365TC124	74
AS828B	52, 75, 78	AS89330	23	IVF1210SS	43	VS365TP124	74
AS828BM	52, 75, 78	AS89800	25	IVF2020	43	VS365TTP124	74
AS828M	52, 75, 78	AS89814	24	IVF2020SS	43		
AS829B-01	52, 75, 78	AS89814H	24	LBE306N	37		
AS830	53, 76, 79	AS89814S	24	LBE307N	37		
AS830B	53, 76, 79	AS89816	24	LBE308N	37		
AS831	53, 76, 79	AS89816H	24	LBE309N	37		
AS831B	53, 76, 79	AS89816S	24	MB4V01	18		
AS832	53, 76, 79	AS89818	24	MBSS001	18		
AS832B	53, 76, 79	AS89818H	24	MR4001	18		



Advancing Healthcare Worldwide®

**Amsino Medical Group Co Ltd.**

708 Corporate Center Drive  
Pomona, CA 91768  
USA

P: +1 800.632.6746

F: +1 909.626.3888

info@amsino.com

[www.amsino.com](http://www.amsino.com)

---

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required. CE mark in the catalog is only based on latest registration status. The CE status may be changed after the latest revision of the catalog.

---

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO®, ACCU-MEASURE®, AMSafe®, AMSafe-3®, AMShield®, AMSure®, AMSure HG®, EZE-VAC®, iRECEPTAL®, One Mask™, Puggie®, R-Polished®, Receptal®, Safgard™, Sal-Sense™, SteriCadi®, Sure-Lok™, THE-POLE-SYRINGE® and VAC-GARD® are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2024, Amsino Medical Group Co. Ltd. All rights reserved.