

RECEPTAL® *Twist*

Safe. Solid. Sealed.

The RECEPTAL® Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options. Engineered with safety and convenience in mind, the intuitive Twist design is simple to use, just one quick twist and done!

- **Premeasured, preloaded and contact-free solidifier bottle** is built into the liner lid for hassle free solidifier administration and waste fluid containment.
- **Integrated solidifier** eliminates the need for open ports, creating a closed system that helps minimize contact with solidifier powder and waste fluid contaminates for better caregiver safety.
- **All-in-one solidifier and liner solution** reduces costs associated with extra inventory needed to carry separate liners and solidifying agents, clean up and potential waste due to solidifier over-measuring.
- **100% leak tested** seamless liners minimize risk of contamination.
- **Latex-free, BPA-free and Non-sterile.**



Maximize safety and minimize costs with a Twist!




Twist solidifier bottle clockwise to release solidifier



After a procedure, simply twist the integrated and preloaded solidifier bottle at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.

Do the Twist with New RECEPTAL® Twist Liners



Ordering Information

Item Code	Description	Qty/Unit
43110S	1000 mL RECEPTAL Twist Liner with VAC-GARD and pre-attached solidifier bottle	40/Case
43115S	1500 mL RECEPTAL Twist Liner with VAC-GARD and pre-attached solidifier bottle	
43116S	1500 mL RECEPTAL Twist Liner with pre-attached solidifier bottle	
43120S	2000 mL RECEPTAL Twist Liner with VAC-GARD and pre-attached solidifier bottle	
43121S	2000 mL RECEPTAL Twist Liner with pre-attached solidifier	

Compatible with RECEPTAL Canisters and Trolley

Item Code	Description	Qty/Unit
43423-01	1500 mL RECEPTAL reusable canister	10/Case
43445-01	2000 mL RECEPTAL reusable canister	10/Case
43476-P2	Adjustable height, Built-in Regulator Manifold Trolley w/ (4) Polycarbonate Canister supports	1/Case



AMSINO

Advancing Healthcare Worldwide®

Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA

+1 800.632.6746 • fax: +1 909.626.3888 • info@amsino.com

For the full RECEPTAL product line, visit us at [Amsino.com](https://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. RECEPTAL and VAC-GARD are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2025, Amsino Medical Group Co. Ltd. All rights reserved.