

AMSure®

Enteral Syringes & Kits with ENFit® Tip



ENFit® Tip

AMSure® Enteral Syringes are compliant with the ENFit® ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections.

AMSure® Thumb Controlled Piston Syringe with ENFit® Tip

AMSure® Thumb Controlled Piston Syringe with ENFit® tip is designed for maximum control allowing for a one-handed administration.



AMSure® Flat Top Piston Syringe with ENFit® Tip

AMSure® Flat Top Piston Syringe with ENFit® tip has visible dosing marks designed for ease of use. Popular choice to deliver viscous fluid with ease.










AMSure® Enteral Feeding Irrigation Kits with ENFit® Syringe

AMSure® Enteral Feeding Irrigation Kit is designed for simplicity and efficiency of care. The non-sterile kit includes a 60 mL thumb controlled piston syringe with ENFit® tip, patient label, 500 mL graduated container. Packaged in a poly pouch or resealable bag.



Order Information

Item Code	Image	Description	Features	Qty/Unit
ENS015		AMSure® 60 mL Thumb Controlled Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS015NS		AMSure® 60 mL Thumb Controlled Piston Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENS115T		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip, transition adapter included	Non-Sterile Resealable Pole Bag*	30/Case
ENS115		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Sterile Poly Pouch	50/Case
ENS115NS		AMSure® 60 mL Flat Top Piston Syringe with ENFit® Tip	Non-Sterile Resealable Pole Bag*	30/Case
ENK126		AMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Controlled Piston Syringe with ENFit® Tip, Patient Label	Sterile Poly Pouch	30/Case
ENK127		AMSure® Enteral Feeding Irrigation Kit Includes: 500 mL Graduated Container, 60 mL Thumb Controlled Piston Syringe with ENFit® Tip, Patient Label	Non-Sterile Resealable Pole Bag*	30/Case

*Resealable pole bag designed for safe and convenient storage.

All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

To view our complete product line,
visit us at www.amsino.com



Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768
800.632.6746 • fax: 909.626.3888 • info@amsino.com

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure® is a registered trademark of Amsino Medical Group Co. Ltd. © Copyright 2024, Amsino Medical Group Co. Ltd. All rights reserved.