

# AMSure®

## Suction Catheters, Trays and Kits

### Suction Catheters

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle tip port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS361

### Suction Catheters and Trays

AMSure Suction Catheter Kits provide hassle-free and convenient extraction of mucus or saliva from the upper respiratory airway.

- **Uniquely designed two-eyed whistle tip offers safe and effective suctioning**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS371

### Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



SCT10

## AMSure® Suction Catheters

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS361</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS361C</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS362</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS362C</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS363</b> ©	10 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS363C</b> ©	10 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS364</b> ©	12 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS364C</b> ©	12 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS365</b> ©	14 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS365C</b> ©	14 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS366</b> ©	16 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS366C</b> ©	16 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS367</b> ©	18 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS367C</b> ©	18 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case

## AMSure Mini Suction Catheter Trays

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>SCT08</b>	Tray, 8 Fr, 17 ½ in. Suction Catheter and Vinyl Gloves	DeLee Tip	100/Case
<b>SCT10</b>	Tray, 10 Fr, 23 in. Suction Catheter and Vinyl Gloves	DeLee Tip	100/Case
<b>SCT12</b>	Tray, 12 Fr, 23 in. Suction Catheter and Vinyl Gloves	Whistle Tip	100/Case
<b>SCT14</b>	Tray, 14 Fr, 23 in. Suction Catheter and Vinyl Gloves	Whistle Tip	100/Case
<b>SCT16</b>	Tray, 16 Fr, 23 in. Suction Catheter and Vinyl Gloves	Whistle Tip	100/Case

## AMSure Suction Catheter Kits

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
<b>AS371</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	1 Vinyl Glove	50/Case
<b>AS381</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	2 Vinyl Gloves	50/Case
<b>AS372</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	1 Vinyl Glove	50/Case
<b>AS382</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	2 Vinyl Gloves	50/Case
<b>AS373</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	1 Vinyl Glove	50/Case
<b>AS383</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	2 Vinyl Gloves	50/Case
<b>AS374</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	1 Vinyl Glove	50/Case
<b>AS384</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	2 Vinyl Gloves	50/Case
<b>AS375</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	1 Vinyl Glove	50/Case
<b>AS385</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	2 Vinyl Gloves	50/Case
<b>AS376</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	1 Vinyl Glove	50/Case
<b>AS386</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	2 Vinyl Gloves	50/Case
<b>AS387</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	18 Fr	2 Vinyl Gloves	50/Case



### Amsino Medical Group Co. Ltd.

708 Corporate Center Drive, Pomona, CA 91768 USA

+1 800.632.6746 • fax: +1 909.626.3888 • info@amsino.com

To view our complete product line, visit us at [www.amsino.com](http://www.amsino.com)

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Items listed with a CE mark are legally available within the European Union (EU) and other regions where CE marking may be required. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSure is a registered trademark of Amsino Medical Group Co. Ltd. © Copyright 2024, Amsino Medical Group Co. Ltd. All rights reserved.